Handi Lift

auto products

WHEELCHAIR LIFT

Sollevatori per disabile

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
2BB 360N

Wheelchair lift

Original instruction translation

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This document is a translation of the manufacturer's original Italian instructions. If further assistance is required, please refer to the original document or forward a request to **Handi Mobil**.

This handbook cannot be reproduce, stored by a filing system or transmitted to third parties, in whole or in part, in whatsoever form or medium, without **Handi Mobil** prior written permission.

Handi Mobil reserve the right to carry out products changes and improvements in order to enhance their quality, at any time and without notice, even if such modifications will not be reported in this publication.

ITALIAN AND FOREIGN PATENT PENDING

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SAFETY RULES AND SYMBOLS

Some precaution are listed in this manual in order to identify potentially dangerous situation in advance and the relevant instructions to avoid them, identified by the following symbols:

 WARNING	This symbol identifies the presence of instructions which need to be read and followed carefully in order to avoid potentially dangerous situations.
---	---

 DANGER	This symbol identifies the presence of essential information needed to avoid potentially dangerous situations that could cause physical injuries and/or damages to the equipment.
--	--

To ensure smooth and safe operation, it is necessary to know the exact procedure for the installation and service of the wheelchair lift.

WARNING

Carefully follow the instructions for the installation and service of the wheelchair lift within this manual. If the instructions within this manual have not been fully understood or further information is required, please contact Handi Mobil immediately. These indications non-fulfilment can cause serious damages, to the equipment and to the vehicles as well as serious injuries to people.

DANGER

Should it be necessary to raise the vehicle being equipped, please check that the characteristics of the lifting device or of the mobile bridge are compatible with the volume and weight of the actual vehicle otherwise there is a risk of serious damage to equipment and vehicles as well as serious injuries to people.

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1. PREFACE

1.1. Classification and general considerations

The wheelchair lift “2BB-360N”, manufactured by **Mobility Networks Ltd**, which is the object of this instruction manual, is in compliance with the safety requirements of 2006/42/CE directive. The “2BB-360N” wheelchair lift in keeping with the aforementioned directive, is to be considered as a ‘machine’ and therefore a declaration of conformity has been released, according to annex II of the 2006/42/CE directive.

All the necessary instructions for the correct installation of the wheelchair lift are detailed for the installer together with a list of maintenance instructions for servicing procedures.

 **WARNING** All installation and service procedures should be carried out in accordance with current health and safety laws.

1.2. Fitter duties / service

Fitter / servicer must:

- carry out the installation in accordance to the indications detailed within this manual;
- carry out the installation by following the indications of the vehicle manufacturer;
- fill and sign the installation conformity declaration;
- carry out the verification check procedure for the first commissioning by following the instructions detailed in this manual;
- register the verification check procedure for the first commissioning in the “use and maintenance” manual;
- give to the final customer all the required documentation by the rules and regulations in force;
- carry out the inspections in compliance with the instructions detailed in this manual;
- register the inspections in the “use and maintenance” manual.

1.3. Documentation

The user should receive with the “2BB-360N” wheelchair lift:

- the original copy of the manufacturer's declaration of conformity;
- the installation conformity declaration both filled and signed by the installer;
- use and maintenance manual filled by the installer.

Handi Mobil is at your complete disposal for further clarification and instructions.

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1.4. Warranty

Handi Mobil assures the hydraulic lift for 24 months from the date of delivery. The warranty covers for defects concerning material quality and product manufacturing.

The warranty does not cover consumables and defects or failures as a result of installation procedure incorrectly carried out or not in compliance with the manufacturer's instructions.

The warranty becomes void in case of impacts caused by accidents and/or tampering carried out by personnel not authorized by **Handi Mobil**.

Handi Mobil refuses all responsibility for damages caused by:

- improper use of the hydraulic lift;
- platform overloading (more than 360kg);
- platform overloading (more than 400kg only FP80150);
- a failure in carrying out "use and maintenance" manual instructions;
- a failure in carrying out maintenance operation as detailed in the "use and maintenance" manual;
- interventions or modifications to the lift without **Handi Mobil** authorization;
- failure in cutting out the electrical circuit, by means of a vehicle battery isolator system;
- lift safety devices removed or tampering;
- wrong manouvres when moving the hydraulic lift;
- presence of objects or people in lift working area.

WARNING

Manufacture's regular inspections dates non-fulfilment will automatically make void the warranty.

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2. GENERAL INDICATIONS

2.1. Technical description

The wheelchair lift "2BB-360N" (Fig.1) is installed on the deck of vehicles used for transporting disabled persons in wheelchairs, allowing them to get in/off the vehicle.

Lift consists of a base fixed to vehicle loading deck a pair of lifting articulation, installed on the sides of the base, and a loading platform, hinged between the articulations.

Deploying/stowing and lifting/lowering movements of the lift are made by means of a parallelogram leverage mechanism driven by a pair of hydraulic cylinders (one for each articulation). The machine is equipped with a hydraulic control unit and an electronic control box which, by means of a remote control, plays the various functional movements; the entire system is electrically powered by the batteries of the vehicle on which it is installed.


Fig. 1

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2.2. Models and installation

The wheelchair lift **"2BB-360N"** is available in three models, which differ from one another in type of loading platform:

- **2BB FP**, with a horizontally folded one platform is deployed automatically during the opening phase and folds automatically during the stowing one, in order to reduce the overall height of the lift at rest, thereby improving the visibility of the driver of the vehicle.
- **2BB P**, with a whole loading platform: during stowing phase platform simply rotates from loading position to vertical one and backwards during deploying phase.
- **2BB SP**, with a loading platform split longitudinally: the two parts of the platforms rotate to form a single loading platform during the deploing phase and they are separated automatically during the stowing one, creating a space that can be used as an emergency exit from the vehicle.

It is possible to prepare the wheelchair lift with following options:

- handrails standar or UK,
- power supply voltage 12V or 24V,
- electro hydraulic power unit right or left.

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
2.4. Overall dimensions

Find following main dimensions of the versions and set-ups available of the wheelchair lift "2BB-360N":

	A	B	C	D	E	F
2BB FP	1106	410	1275	1233	800	1500
	1106	410	1275	1233	840	1380
	1006	410	1175	1233	740	1380
	966	355	1115	1074	700	1150
2BB P	1106	305	1275	1480	900	1380
	1006	305	1175	1402	800	1300
	966	305	1115	1300	760	1200
2BB SP	1106	500	1275	1600	900	1500
	1106	500	1275	1480	900	1380
	1006	450	1175	1400	800	1300
	966	430	1115	1300	760	1200

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2.5. Safety devices

The rules described above already provide a good level of safety, capable of avoiding danger if followed properly.

The “**2BB-360N**” lifts are also equipped with additional safety and security devices:

- **Safety handrails:**
ensure the function of “firm hold” during the lifting / lowering phase.
- **Internal flap:**
to let wheelchair get in / off the platform; it is mounted in order to stay vertical when the platform is not at vehicle deck level in order to avoid crushing and shearing risk.
- **External flap:**
to let wheelchair get in / off the platform; it is mounted in order to stay vertical when the platform is not at ground level in order to avoid the risk of wheelchair falling from the platform of the wheelchair.
- **Protection devices (fitted on moving parts):**
they are carter with the function of protection from the moving parts.
- **Protection against overturning:**
stops platform lowering in case of an obstacle which induces overturning.
- **Safety pressure switch:**
stops stowing when platform is still loaded.
- **Protection against overloading:**
blocks lifting in case there is a load greater than 125% of the nominal loading capability of the lift (360kg/400kg).
- **Protection against hydraulic system leakage:**
avoid platform descent speed from exceeding operating one by 50%.
- **Safety hook:**
avoid the unintentional opening of the platform in case of low pressure in the hydraulic system when the lift is closed.
- **Overall dimensions signalling:**
marker strips are stuck on the perimeter of the lift; besides intermittent led lights (optional) are provided, screwed to the top front portion of both the lifting articulation.

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
Emergency maneuvers safety sticker positioned on the front of the lifting articulation on the side of the hydraulic control unit.

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Lift safety use sticker, positioned on the front of the lifting articulation on the opposite side of the hydraulic control unit.

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Identification plate riveted on internal side of lift base.

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
Safety sticker on the internal flap.

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
Identification sticker of the connector of the wired remote control panel placed on the front of the control unit box.

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Identification sticker of the joystick command (if present) for the down/up maneuvers placed in correspondence to the related command.

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

Identification sticker of the manual pump placed on the control unit box corresponding with the slot for the hand pump lever.

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Identification sticker of the manual valve placed on the control unit corresponding to the return manual valve for descent.

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Identification sticker of the connectors for the automatic opening and closure of the doors and of the remote control placed at the bottom and on the front of the control unit box alongside the relevant connectors, if present.

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WARNING In the case of loss or excessive wear, identification tags and safety stickers should be rapidly replaced. Please contact an authorized Handi Mobil centre for details.

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3. PRELIMINARY OPERATIONS

3.1. Transport and movement

For transportation, the lift must be secured to a pallet and packed with sheets of protective plastic and cardboard. For transportation purposes, a fork lift truck or a hoist is advised.

All packing materials are ecological and recyclable. The packing materials should be disposed off correctly. The local council or refuse department will be able to offer advice regarding disposal requirements.

 **DANGER** Cardboard and protective plastic sheets used for packing purposes can cause suffocation.

3.2. Package contents

The size of the wheelchair lift package "2BB-360N" is 135x42 cm.

The height of the wheelchair lift package varies according to the model type.

140 cm for the **2BB FP** model;

170 cm for the **2BB P** and **2BB SP** models.

The package weight is 13 kg.

 **WARNING** During the unpacking process be very careful to not damage the contents.

 **WARNING** Before beginning any installation procedure of the wheelchair lift, you should:

- verify if there are any specifications from the vehicle manufacturers to respect;
- remove from the vehicle any unnecessary object that could impede the installation procedure of the lift (spare wheel, accessories etc);
- disconnect the electrical supply from the battery;
- disconnect any electronic control units from the vehicle as specified by the manufacturer.

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4. STABILITY VERIFICATION

Before starting to position and fix the lift, as described in *Chapter 5*, it is necessary to verify the stability of the vehicle with the lift installed open and maximum overhang.

Here are the instructions for this verification procedure:

- Obtain from the vehicle registration documents, the information required for stage (**P**) and the rear overhang limits (**S**) and write them in *Section A* of *Chapter 1* of "use and maintenance" manual ("Log verification for the first commissioning");
- From the vehicle obtain measurement (**A**) from the loading base to the height of the bumpers (**Fig.2**) and compare it to the values as detailed in the table below, obtain measurement (**B**) from the distance of the external edge of the lift base to the external edge of the bumpers; the values of **B** as detailed in the table are valid for all models of lift apart from 76120 and 70115 dimensions.
- The values for these 2 models are detailed as **B*** in the table; transcribe the values obtained of **A**, **B**, **B*** in *Section A* of *Chapter 1* of "use and maintenance" manual ("Log verification for the first commissioning");


Fig. 2

A	(mm)	0	20	40	60	80	100	120	140	160
B	(mm)	210	230	245	265	280	295	305	320	330
B*	(mm)	230	250	270	290	300	315	315	315	315


Fig. 3

	Platform dimensions [cm]			d1(m)	G1 (daN)	d2(m)	G2 (daN)	M(daN)
	Width	X	Length					
2BB FP	84	X	138	0,10	90	1,12	70	160
	80		150	0,10	90	1,12	74	160
	74	X	138	0,11	88	1,12	62	150
	70	X	115	0,95	80	1,00	52	132
2BB P	90	X	138	0,10	90	1,12	57	147
	80	X	130	0,11	88	1,12	49	137
	76	X	120	0,95	80	1,00	45	125
2BB PS	90	X	138	0,10	90	1,12	70	160
	80	X	130	0,11	88	1,12	62	150
	76	X	120	0,95	80	1,00	52	132

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- By entering the values previously obtained into the following formula you will be able to obtain the minimal reaction (R_{min}) for the anterior axle of the vehicle and transcribe it in *Section A* of *Chapter 1* of "use and maintenance" manual ("Log verification for the first commissioning");
- Obtain from the vehicle registration documents, the value of the ultimate load capacity for the anterior axle (R) and transcribe it in *Section A* of *Chapter 1* of "use and maintenance" manual ("Log verification for the first commissioning");
- If the value of R is higher than the value of R_{min} , the stability of the vehicle is verified; complete *Section A* of *Chapter 1* of "use and maintenance" manual ("Log verification for the first commissioning"), with signature and stamp and proceed with the lift installation procedure as described in the following chapter.

Stability check is positive when: $R > R_{min}$

WARNING Please follow the instructions in this manual carefully when carrying out the stability verification of the "2BB-360N" wheel chair lift. If the instructions in this manual have not been properly understood or if further information is required, please contact Handi Mobil immediately. Serious damage to vehicles and human injury could be caused if these instructions are not carefully fulfilled.

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5. LIFT POSITIONING AND FIXING

Here are all the necessary instructions for the installation of the **"2BB-360N"** wheelchair lift for the vast majority of vehicles on sale. Not being able to cover all possible cases, we advise you to contact **Handi Mobil** for any further information.

- Pressurize the hydraulic system of the lift by using its hand pump as follow: screw the provided lever in its proper place in the hydraulic control unit (**Fig.4b**) and put the hand pump into action with an alternate vertical movement until the pump's resistance blocks the lift in the "closed" position.
- Open the door(s) of the vehicle, where the lift will be placed, and block them on the "open" position". Measure the height and the width of the compartment and verify that they are bigger than the overall dimensions of the lift (**Fig.5**).
- Using a forklift or an equivalent means with a bigger capacity than the total weight of the lift (see technical specifications), place it at the same height as the surface of the vehicle (placing compartment). Push the lift inside the vehicle in a central position in respect to the volume of the compartment. Align the external edge at the base of the lift parallel to the closing edge of the door(s).
- Position the lift in order to its base external edge is at distance "B" (calculated in *Chapter 4*) from vehicle bumper external edge.
- Verify, both on the inside and the outside of the vehicle, that the doors of the vehicle close correctly, with no interference from the lift. Measure the minimum distance between the doors and the lift; if it is more than 40 mm (minimum distance) (**Fig.6**) it is possible to draw up the lift to the doors until it reaches the minimum distance. Pay attention to keep the external edge at the base of the lift and the closing edge in a parallel position. Also pay attention to the centering of the lift inside the width of the compartment.


Fig. 4

Fig. 5


Fig. 6

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WARNING If lift is placed at a distance lower than "B" quote (calculated in Chapter 4) stability verification must be made again and "use and maintenance" manual Chapter 1, section A must be updated ("Log verification for the first commissioning").

- Before fixing the lift, check that there are no interferences (in correspondence to the tightening screws) with the parts under the frame, such as hydraulic conduits or wireways, electrical wiring, cables etc. Move the lift sideways to avoid interference.
- After having identified the exact position of the lift on the fitting surface area of the vehicle, use the eight holes at the base of the lift as a guide for the drill holes to be made on the surface of the vehicle.
- It is important to ensure that the supporting surface to which the lift will be fastened, is stable and that the holes do not affect the stability of the lift. The lift is provided with fastening plates and these have to be positioned, in correspondence to the holes, under the frame (Fig.7). They need to be tightened using n°8 M10 screws with their screw nuts and washers (Fig.8) using a torque equal to 50 Nm.

WARNING If it is essential to use screws with a length greater than those of the screws provided, make sure that they are M10 class 8.8 screws.

WARNING Most vehicles on sale have a corrugated structure (Fig.9) and in these cases the empty spaces need to be filled with plates which act as spacers in the bolt fixing area.


Fig. 7

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
Fig. 8

WARNING If, in order to stiffen the area that will support the lift, welding is needed, first unplug all the existing electrical connections on the vehicle and follow the manufacturer's instructions with care.

WARNING It is the responsibility of the installer to verify the adequate resistance and crush proof characteristics of the surface of the vehicle to which the lift will be attached.

Handi Mobil declines all responsibility for any damage caused by such instructions and verifications non-stilfilment and it will automatically consider void the warranty of the lift.


Fig. 9

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6. ELECTRICAL CONNECTIONS

6.1. Power supply

After a correct positioning and fastening of the lift to the loading deck of the vehicle, connect the eyelet connectors of the 16 mm² wires to the battery of the vehicle, in particular, connect the red wire to the positive pole and the black wire to the negative pole.

 **WARNING** Considering the 12 V or 24V application, connect the earth wires directly to the negative pole of the battery of the vehicle at the motor gear of hydraulic power unit.

 **WARNING** If wires need to go through walls or panels use specific rubber gromits to avoid damaging sheaths and wires.

 **WARNING** Connecting procedures must comply with the instructions given out by the vehicle manufacturers.

 **WARNING:** The connection to the vehicle batteries must be carried out by passing towards the vehicle electric battery-isolator.

Handi Mobil declines all responsibility for any damage caused by electric connections not in compliance with what has been stated in this manual and with the technical specifications of the vehicle manufactures.

6.2. Output signal

The lift is provided with a six pins female connector for signal wires not connected; in this situation the lift is not enabled to function.


Before installation, if any verifications of the lift are necessary, connect to the power supply and connect the six-way female connector to its correspondent male connector which is placed inside the hydraulic control unit.

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During installation, remove the bridge on the female connector and perform the connections as shown in under picture, installing a warning light inside the cabin, a beeper and a system enabling switch; connect the female connector to its male connector positioned inside the hydraulic control unit.

Handi Mobil decline all responsibility for damage caused by non fulfillment of these instructions and automatically deem the warranty to be void in such cases.

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7. HYDRAULIC SYSTEM

7.1. Hydraulic unit

The hydraulic unit of the lift can be delivered in two versions: manual controls on the left side (left version, see **Fig. 10**) or manual controls on the right side (right version); both versions are suitable only for vertical installation.

- 1) Electro-pump electric engine
- 2) Seat of manual pump
- 3) Adjustable throttle
- 4) Knob lowering electro valve
- 5) Contactor
- 6) Manual pump
- 7) Cap
- 8) Oil tank


Fig. 11

In order to assure to the electro-pump "1" (**Fig. 11**) an optimally flow of hydraulic oil, it is necessary that the level oil in the tank "8" (**Fig. 11**) it is, with platform stowed, above of the minimum level.

It is possible to check the hydraulic oil level into the tank by unscrewing and removing the cap "7" (**Fig. 11**) with the platform stowed.

⚠ WARNING When the vehicle engine is switched off, it is suggested not to operate the hydraulic unit for more than one minute at a time, in order not to overload batteries.

⚠ WARNING When filling up oil, do not mix oils having different features and gradations.

⚠ WARNING When checking and filling up oil, keep lift stowed.

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7.2. Pressure socket

WARNING

These instructions are for the exclusive use of appropriately trained technical personnel.

DANGER

Before using the pressure socket make sure that platform is firmly on the ground in order to get the system pressure at zero.

Inside the hydraulic control unit there is a pressure socket. To measure the pressure in the hydraulic circuit; use a 6 allen wrench to screw off the lid of the pressure socket placed on the T-piece found on the top of the distribution block of the control unit, in correspondence to the pressure switch and connect through the male connector 1/4" .

WARNING

After disconnecting from the pressure socket, make sure to screw and tighten the lid.

Handi Mobil decline all responsibility for damage caused by non fulfillment of these instructions and automatically deem the warranty to be void in such cases.

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8. VERIFICATION FOR FIRST COMMISSIONING

The following instructions are necessary for the first commissioning of the “**2BB-360N**” wheelchair lift after it has been installed on the vehicle. Note that every record is equipped with a description of the verification check, time estimations, number of personnel required and necessary materials and equipment.

After each verification has been carried out, complete, sign and stamp *Section B* of *Chapter 1* of “use and maintenance” manual (“Log verification for the first commissioning”).

WARNING

These instructions are for the exclusive use of appropriately trained technical personnel.

WARNING

For the verification check of first commissioning of the “**2BB-360N**” wheelchair lift the instructions detailed in this manual must be carefully fulfilled.

If the instructions in this manual have not been properly understood or if further information is required, please contact **Handi Mobil** immediately. Serious damage to vehicles and human injury could be caused if these instructions are not fulfilled.

Handi Mobil declines all responsibility if the verification check for the first commissioning is not followed properly and if the appropriate records are not compiled correctly (in *Section B* of *Chapter 1* of “use and maintenance” manual) and will automatically consider the warranty void.

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9. STANDARD INSPECTIONS

Standard inspections should be carried out 6 month by **Handi Mobil** Authorized Service Centre. They allow for greater operational safety and guarantee a longer lifetime for the **2BB-360N** wheelchair lift.

Here are the instructions for standard lift inspections to be followed by service centre. In particular, in every inspection form you can find a description of the verification carried out, and of the possible intervention in case of a negative check test result, time estimations, number of personnel required and necessary materials and equipment.

The record for the standard inspection, present "use and maintenance" manual, must be filled, stamped and signed after every test check.

 **WARNING** These instructions are for the exclusive use of appropriately trained technical personnel.

 **WARNING** In the standard inspection phase, the servicer is also required to carry out any maintenance interventions even if they are the responsibility of the user of the lift, and as described in *Chapter.7*.

 **WARNING** After **2000 HOURS** of functioning and always after **10 YEARS**, the lift should undergo a full and complete service at a **Handi Mobil** authorized centre, who will subsequently issue a new user manual with new inspection procedures prior to resuming service.

IF THIS IS NOT THE CASE THE LIFT MUST NOT BE OPERATED.

 **WARNING** For the standard inspections on behalf of the servicer of the **2BB-360N** wheelchair lift the instructions detailed in this manual must be carefully adhered to.

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
Find following “Handi Lift” operation system described in consideration of micros M1, M2 ed M3 positions:

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10. ADDITIONAL INSPECTIONS

The authorized **Handi Mobil** , service centre which has been contacted by a customer encountering anomalies or malfunctioning of the wheelchair lift, must carry out the required checks and if necessary operates in order to solve the associated problems. At the end of the additional inspection service centre must fill, stamp and sign the relevant form in lift “use and maintenance” manual.

WARNING

In the additional inspection phase, the service centre is also required to carry out any maintenance interventions even if they are the responsibility of the user of the lift, and as described in *Chap. 7*.

WARNING

For the additional inspections on behalf of the servicer of the **2BB-360N** wheelchair lift the instructions detailed in this manual must be carefully adhered to.

If the instructions in this manual have not been properly understood or if further information is required, please contact **Handi Mobil** immediately. Serious damage to vehicles and human injury could be caused if these instructions are not adhered to.

Handi Mobil decline all responsibility for damage caused due to incorrect maintenance or lack of maintenance of the lift and automatically deem the warranty to be void in such cases.

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Fig. 14

11.1. KIT LIGHTS

The Light kit is composed of n°2 lintermittant led lights, n°4 tapping screws and n°2 electrical connectors.

To install the Light kit follow the instructions (on each side of the lift).

- Thread the led light power cord into the hole (Fig. 12).
- Fix the led lights to lift by means of the n°2 tapping screws (Fig. 13).
- Connect the lights to the branch of electrical system, which is identified with the label, 'lights' (Fig. 14 - Fig. 15).

Fig. 13

Fig. 15

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
Fig. 16

11.2. KOYSTICK CONTROL KIT (SEC0008)

The Joystick control kit is composed of n°1 wired Joystick divided into 3 parts and with cap, n°1 cable conduit 250 mm in length, n°1 cable gland and conduit 600 mm in length, n°1 4 pins male connector and n°1 4 pins female connector, n°4 female Faston connectors, n°4 male Faston connectors and n°1 sticker.

To install the Joystick control kit on the 2BB-360N, follow the instructions below:

- Remove the cap from the left side hand rail.
- Fit the 250 mm conduit inside the left vertical arm and the 600 mm conduit inside the lower left horizontal arm.
- Feed the Joystick side cable through the hand rail, without re positioning the cap, feed the cable as shown in Fig. 14 inserting into the 250 mm conduit.
- Feed the central cable into the 600 mm conduit.
- Feed the free end under the base until it reaches the "JOY" cable of the control unit.
- Feed the free end under the base until it reaches the control unit.
- Connect the electrics as detailed in diagram (Fig. 16).
- Test the switch and rotate the cap until the command position for 'up' is towards the top.
- Apply the cap and the sticker (Fig. 17).
- Check for the correct positioning of the cable when the lift is in operation.

Fig. 17

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TO REFER A SEC0008 WIRING (COMMERCIAL MODEL)

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Fig. 18

11.3. JOYSTICK CONTROL KIT (SEC0007)

Joystick control kit is composed of n°1 wired joystick and cap, n°1 cable conduit, n°1 4 pins connector and n°1 sticker. To install the Joystick control kit on Handi Lift (Mobility model), follow the instructions below:

- Remove the cap from the left side hand rail.
- Fit the cable conduit as shown in Fig.18.
- Feed the cable into the handrail and feed through, without re positioning the cap as shown in Fig.19 and feed it into the conduit.

Fig. 19

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- Fit the 4 pins connector (from the pack) to the free end of the cable and connect it to the branch of the lift's electrical system identified with the label 'Interr'.
- Test the switch and rotate the cap until the command position for 'up' is towards the top.
- Apply the cap and the sticker (Fig.20-21).

Fig. 20

Fig. 21

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Fig. 24

12.1. LONGITUDINAL CHECK

Set to ZERO the measuring instrument and put it on the work bench (or vehicle deck for the installer) where the lift base is fixed.

1. Raise platform to vehicle deck level (Fig. 22);
2. Place the measuring instrument on the platform in the longitudinal direction; then check that angle measurement is between $+ 1.5^\circ$ and $+1.7^\circ$ (Fig. 23 e Fig.3).

Fig. 23

3. If measurement is out of the range of point 2, adjust platform inclination by means of platform adjusting screws (Fig. 25) :
4. Loosen the screws (first right one and then the left one) by means of a measure 6 allen spanner (fig. 4); turn them clockwise to get a positive angle or counterclockwise to get a negative angle.
5. Do again the check and verify that the new angle is in point 3 range.


Fig. 25

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Fig.3

3. Loosen the chains by means of the 2 setting nuts one turn unscrewing (Fig.3).
4. Unscrew the setting tilt platform headless screws (Fig.4 and after paragraph) until two platform sides are realigned.
5. Check external flap correct opening at the ground level; if it is rightly then reassemble all components, else restart procedure from point 1.


Fig.4

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
12.3. SAFETY PRESSURE SWITCH ADJUSTMENT

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1. Since platform at ground level, put on its edge 25 kg payload (see picture 1).
2. Play "Up/Stow" button to raise platform.
3. If platform (with 25 kg payload) doesn't stop itself at vehicle floor level (see picture 2), starting closing, and hydraulic pump stays on, pressure switch adjustment is required in compliance with instructions starting from step (4); on the contrary go to step (6).
4. Unfit hydraulic unit cover, by unscrewing its n.5 screws,
5. Insert in pressure switch socket measure 1,5 (2 for old models) allen spanner and, while playing "Up/Stow" button, turn it a little bit counterclockwise (see picture 3) until platform stops moving and hydraulic pump switches off.
6. Remove the payload and check that, playing "Up/Stow" button, platform keeps on closing until it reaches final stowing position (see picture 4). If check is positive, go to step (9), on the contrary go to step (7).

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7. If not already done, unfit hydraulic unit cover, by unscrewing its n.5 screws,
8. If not already done, insert in pressure switch socket measure 1,5 (2 for old models) allen spanner and, while playing "Up/Stow" button, turn it a little bit clockwise (see picture 3) until platform is completely stowed.
9. Once adjusted, repeat movements for five times to check the correct setting (platform should stop all the time at vehicle floor level, if not, repeat this procedure starting from step 3).

At the end of adjustment, check that allen spanner is been unfitted and hydraulic unit cover is properly mountend by means of its five screws.

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13. TROUBLESHOOTING


Find following troubleshooting table for all the failure not listed or caused by mechanical or hydraulic damages, please, contact as soon as possible the nearest service centre.


Trouble	Possible Cause	Solutions
1 Cabin warning light OFF 	System enabling switch OFF	 Turn ON system enabling switch
	Hydraulic power unit six-poles connector not plug	 Plug six-poles connector
	Manual reset fuse open	 Reset the fuse
	Main power fuse blown	 Replace main power fuse
	cabin warning light broken	 Replace cabin warning light
	electronic control unit fuse blown	 Replace electronic control unit fuse
2 Platform does not deploy 	Safety hook in locking position	 Stow the platform
	Wired remote control not properly connected	 Connect properly wired remote control
	Lowering valve solenoid electrical circuit open	 Close solenoid electrical circuit
	Lowering valve solenoid broken	 Replace solenoid lowering valve
	Lowering valve broken	 Replace lowering valve
	Wired remote control broken	 Replace wired remote control
	Electronic control unit broken	 Replace electronic control unit
3 Platform does not move aligned 	Throttle valve broken (one of the two)	 Replace both throttle valves
	Lifting cylinder sealings damaged	 Replace lifting cylinder sealings
	Lifting cylinder rod bended	 Replace lifting cylinder
4 Platform does not stop at loading deck during deploying 	Micro-switch M1 not properly connected	 Connect properly micro-switch M1
	Micro-switch M1 electrical circuit open	 Close micro-switch M1 electrical circuit
	Micro-switch M1 cam not properly adjusted	 Adjust properly micro-switch M1 cam
	Micro-switch M1 broken	 Replace micro-switch M1
5 Platform does not make lowering movement 	Wired remote control not properly connected	 Connect properly wired remote control
	Micro-switch M2 not properly connected	 Connect properly micro-switch M2
	Micro-switch M2 checking plate not properly fixed	 Fix properly micro-switch M2 checking plate
	Micro-switch M2 electrical circuit open	 Close micro-switch M2 electrical circuit
	Lowering valve solenoid electrical circuit open	 Close solenoid electrical circuit
	Micro-switch M2 broken	 Replace micro-switch M2
	Lowering valve solenoid broken	 Replace solenoid lowering valve
	Lowering valve broken	 Replace lowering valve
	Wired remote control broken	 Replace wired remote control
	Electronic control unit broken	 Replace electronic control unit
6 Platform goes down in uncontrolled way 	Lowering valve knob not tighten	 Tighten lowering valve knob
	Lowering valve broken	 Replace lowering valve
7 Platform deploys/goes down too fast 	Adjustable throttle valve not properly adjusted	 Adjust properly adjustable throttle valve
	Adjustable throttle valve broken	 Replace adjustable throttle valve
8 Platform does not raise 	Lowering valve knob not tighten	 Tighten lowering valve knob
	Lowering valve broken	 Replace lowering valve
	Hydraulic pump teleruptor broken	 Replace hydraulic pump teleruptor

Trouble		Possible Cause	Solutions
		Hydraulic pump damaged	 Replace hydraulic pump
		Wired remote control broken	 Replace wired remote control
		Electronic control unit broken	 Replace electronic control unit
9 Platform does not stop at loading deck at the end off lifting		Micro-switch M3 not properly connected	 Connect properly micro-switch M3
		Micro-switch M3 electrical circuit open	 Close micro-switch M3 electrical circuit
		Micro-switch M3 cam not properly adjusted	 Adjust properly micro-switch M3 cam
		Micro-switch M3 broken	 Replace micro-switch M3
10 Platform does not stow		Lowering valve knob not tighten	 Tight lowering valve knob
		Pressure switch not properly set up	 Set up properly pressure switch
		Lowering valve broken	 Replace lowering valve
		Hydraulic pump teleruptor broken	 Replace hydraulic pump teleruptor
		Hydraulic pump damaged	 Replace hydraulic pump
		Wired remote control broken	 Replace wired remote control
		Electronic control unit broken	 Replace electronic control unit
11 Platform does not stow entirely		Pressure switch not properly set up	 Set up properly pressure switch
		Electronic control unit delay parameter uncorrect	 Change electronic control unit delay parameter

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
By:

Mobility Networks Holdings Ltd

**12 Estuary View Business Park, Whitstable, Kent,
CT5 3SE, UK**

**Tel: +44 (0)1227 505022 - Fax: +44 (0)1227 505023 -
www.mobilitynetworksgroup.com**

**Mobility Networks is a trading name of Mobility
Networks (Holdings) Ltd.**

Registered in England & Wales number 08295602